

Group 3 PET Tool Script

URL: <http://g3cap.tk/staging>

Admin Credentials:

Email: g3cap.tk@gmail.com

Pass: Capstone1 (Case Sensitive)

Script:

Normal User:

1. Go to <http://g3cap.tk/staging>
2. Register an account:
 - a. Fill out form and Submit
 - b. Check email for account validation
3. Return to login page
4. Sign In with new account
5. Search Known Phages:
 - a. Select multiple phages/clusters/subclusters/enzymes and click search
 - b. Repeat Step 5 until you are satisfied with the results
6. Unknown Phage Search:
 - a. Click the checkbox besides the "Search" button that reads "Add Unknown Phage?"
 - b. Fill out the corresponding form with unknown phage data.
 - c. Click "Search"
 - d. Compare the phage data to the unknown
7. Licensing Information:
 - a. Click the License Page in the Navbar at the top.
 - b. Follow the link to the Licensing Information.
8. Demographics:
 - a. Click the Demographics Page in the Navbar at the top.
 - b. Choose pie graph data to see different visualizations.
9. Log Out by clicking the Log Out Button at the top right in the Navbar.

Admin:

1. Sign in with admin credentials
2. Click on the Admin Tools Page in the Navbar at the top.
3. Adding a phage:
 - a. Enter text for each into the text fields
 - b. Click "Submit" on the right
4. Deleting an enzyme:
 - a. Select enzyme for deletion from the select box

- b. Click "Delete Enzyme"
- 5. Delete an account:
 - a. Click "Account Management" Tab near the top
 - b. Click the account to delete
 - i. (CAREFUL!! THERE IS CURRENTLY NO SAFEGUARD AGAINST ACCOUNT DELETION. SELECTING AN ACCOUNT ON THE TABLE WILL IMMEDIATELY DELETE THE ACCOUNT)
 - ii. If you need to add another admin account, currently it has to be done by a database administrator.

Visit <http://g3cap.tk/staging/README> for helpful documentation.